ChromoLympho-B MIX

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Instruction for use

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Reagent for the *in vitro* stimulation of human B-lymphocyte.**Material supplied**

Product	Code	Format	Storage
ChromoLympho-B MIX	EKAMP010M	50 test/50 tubes	-18/-22°C

Each tube contains ChromoLympho-B Proliferation Factor (synthetic oligodeoxynucleotide, 27bases) and IL-2 Human (Recombinant Human Interleukin-2, source: Yeast) premixed and lyophilised.

Storage and Validity

ChromoLympho-B MIX is shipped at room temperature and must be stored at -18/-22°C upon arrival. Lyophilized powder is stable at -18/-22°C until expiration date provided on the tube. After reconstitution, store liquid working aliquots at -18/-22°C until use.

ChromoLympho-B MIX solution is stable during 6 months at -18/-22°C upon reconstitution in sterile conditions.

Please avoid freeze-thaw cycles.

NOTE: aspect and amount of the white pellet in the tube may vary from tube to tube, as a result of salts precipitations during lyophilization step.

General protocol (for reference only)

In vitro stimulation of human B-lymphocytes, in cultured peripheral blood cells or cultured bone marrow cells from patients with chronic lymphocytic leukaemia (CLL) and other B-cell lymphoproliferative syndromes.

1. Prepare cells in 9 ml of culture medium ($1-2 \times 10^6$ cells/ml; blood or bone marrow).
2. Before opening the tube of ChromoLympho-B MIX, do a short spin of 20 seconds to collect the pellet down.
3. Dissolve lyophilized powder in 1 ml of culture medium (Chromosome M*) and ensure solution is uniform by repeated pipette mixing.
4. Add 1 ml of MIX into 9 ml of culture cell preparation.

* Chromosome M: complete culture medium for bone marrow samples, sold by Euroclone S.p.A.

General protocol for Synchronized cultures (for reference only)

1. After 48-72 H of culture, add 0.2 ml of Synchroset solution A (Ref. EKAMTS008**) into the 10 ml cell culture.
2. After overnight incubation, add 0.2 ml of Synchroset solution B (Ref. EKAMTS008), without changing the culture medium, and incubate further 5 hours.
3. Add few drops of Colcemid ($10 \mu\text{g/ml}$ - Ref. EK0041B***), and incubate at $37 \pm 0.5^\circ\text{C}$ during 60 minutes (if Synchroset is used for the first time, different incubation times with Colcemid should be tried, i.e. 30 min, 45 min, 1h).
4. Follow the standard protocol for metaphase preparation.

****EKAMTS-008:** Synchronization kit sold by Euroclone S.p.A.

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*Chromosome M: terreno di coltura completo per campioni di Midollo Osseo, commercializzato da Euroclone S.p.A.
**EKAMTS008: Synchronization kit, commercializzato da Euroclone S.p.A.
***EK0041B: Colcemid commercializzato da Euroclone S.p.A.
*****EK0041B:** Colcemid sold by Euroclone S.p.A.