

## FISH Reagent Kit

**REF** KBI-60005

Containing sufficient volume for up to 25 or 5x5 tests

Kreatech Biotechnology B.V.  
Vlierweg 20  
1032 LG Amsterdam  
The Netherlands  
[www.LeicaBiosystems.com](http://www.LeicaBiosystems.com)

### 1 Intended use/purpose

Fluorescence in situ hybridization (FISH) identifies or labels target genomic sequences so that their location can be studied. DNA sequences from appropriate chromosome specific probes are first labeled with reporter molecules. The labeled DNA probe is then hybridized to the metaphase chromosomes or interphase nuclei on a slide. After washing, the specimen is screened for the reporter molecules by fluorescence microscopy.

For use on metaphase and interphase cells from peripheral blood cultures or direct preparations prepared by standard cytogenetic methods.

This pretreatment kit is specifically developed to obtain optimal results on FRESHLY prepared cytological samples.

This product is intended for In Vitro Diagnostic Use.

### 2 Product Description

This pretreatment kit is for non-automated FISH hybridizations specifically developed to use for pretreatment and enzyme digestion of metaphase and interphase cells from peripheral blood cultures or direct preparations prepared by standard cytogenetic methods before DNA FISH hybridization procedure. This kit is specifically developed to obtain optimal results on FRESHLY prepared cytological samples. The use of the kit will reverse the cross-links in the cells, permeabilize the cell membrane and digest cell proteins to ease the binding of the FISH probes with the specific cell DNA.

### 3 Instructions for use

#### 3.1 Specimens

For use on metaphase and interphase cells from peripheral blood cultures or direct preparations prepared by standard cytogenetic methods.

This pretreatment kit is specifically developed to obtain optimal results on FRESHLY prepared cytological samples.

Specimens pre-analytics: To obtain a successful FISH the pre-analytics play a major role. The specimen should be correctly fixated, aged and pretreated according to a well-established protocol. Over-fixation, under fixation or poor pretreatment of specimen can cause failure of the FISH procedure or at least low reproducibility/signal intensity. Kreatech FISH probes have been tested on metaphase/interphase nuclei derived from cultured peripheral blood and/or bone marrow. Slides should be stored according to the recommended storage and handling conditions for optimal results, see Leica website.

#### 3.2 Pretreatment

Protocol below is for the use of 5 slides in an upright coplin jar of 50 ml.

1. Pretreat freshly prepared sample slides in 50 ml of 2 x SSC / 0.5% Igepal (LK-105B) at 37 °C for 15 minutes (min).
2. Dehydrate in 70%, 85% and 100% ethanol for 1 minute each. Airdry. Proceed with **Probe Preparation**.

#### 3.3 Probe preparation:

Kreatech™ FISH probes are supplied Ready to Use (RtU). SE, ST, and WC Kreatech™ FISH probes are provided 5 x concentrated and must be diluted as follows: 2 µl 5 x conc. Probe in 8 µl FISH Hybridization Buffer (FHB or WCB, supplied with probes). To combine several 5 x conc. probes, replace FISH Hybridization Buffer (FHB or WCB) with 2 µl for each probe added.

### 3.4 Co-denaturation

Apply 10 µl of probe per 22 x 22 mm field. Cover with glass cover slip and seal with Fixogum or rubber cement. Denature sample and probe on a ThermoBrite (TS-01/02) at 75 ±1 °C for 5-10 min. Continue with **Hybridization**.

### 3.5 Hybridization:

Incubate overnight at 37 ±1 °C in a ThermoBrite (TS-01/02) or in a humidified chamber.

### 3.6 Post-hybridization Wash:

1. Pre-warm 50 ml Wash Buffer I (0.4 x SSC / 0.3% Igepal) (LK-102) to 72 °C
2. Remove rubber cement.
3. Place slides in 50 ml of Wash Buffer II (2 x SSC / 0.1% Igepal) (LK-103), incubate for 2 min. at RT. Slide off coverclips. Re-use the wash buffer only once for a total of 10 slides.
4. Place slides in 50 ml of pre-warmed Wash Buffer I (0.4 x SSC / 0.3% Igepal) (LK-102), incubate for 2 min at 72 °C (+- 1 °C) without agitation. Re-use the wash buffer only once for a total of 10 slides.
5. Place slides in 50 ml of fresh Wash Buffer II (2 x SSC / 0.1% Igepal) (LK-103), incubate for 1 min. at RT without agitation. Re-use the wash buffer only once for a total of 10 slides.
6. Dehydrate in fresh 70%, 85% and 100% ethanol, incubate for 1 min each at RT. Air dry at RT and proceed to **Counterstaining**.

### 3.7 Counterstaining:

Apply 15 µl DAPI counterstain (DAPI/Antifade 0.1 µg/ml) (LK-095B) or a dilution of this using counterstain diluent (LK097B) and apply glass cover slip. Proceed with microscopy.

### 3.8 Procedural recommendations:

Optimal prepared and pre-treated cytological slides are key to obtain a successful FISH result. Under or over pre-treatment does not only influence the morphology but also the signal strength (and thus the ease of analysis).

Temperature and buffer concentration (stringency) of hybridization and washing are important, as lower stringency can result in non-specific binding of the probe to other sequences, and higher stringency can result in a lack of signal. Incomplete denaturation of target DNA can result in lack of signal.

### 3.9 Recommendations for fluorescence microscopy

For optimal visualization use a well maintained and regularly calibrated microscope equipped with a 100W mercury lamp or other appropriate light source (e.g. Metal Halide Bulb (120W X-cite) or LED based light sources) and a 63x or 100x fluorescence objective. Triple band-pass filters (DAPI/FITC/Cy3 or DAPI/FITC/TRITC) are used to view multiple colors, single band-pass filters are used for individual color visualization. See the Leica website for recommendations on appropriate filters for each microscope brand: <http://www.leicabiosystems.com/ihc-ish/kreatech-fish-probes/>

Suitable excitation and emission range for Kreatech™ fluorophores:

| Fluorophore         | Excitation | Emission   |
|---------------------|------------|------------|
| DAPI                | 360 ±20 nm | 460 ±30 nm |
| PlatinumBright™ 415 | 429 ±15 nm | 470 ±15 nm |
| PlatinumBright™ 495 | 495 ±10 nm | 517 ±10 nm |
| PlatinumBright™ 550 | 550 ±10 nm | 580 ±10 nm |

*The shown bandwidth corresponds to 80% of the maximum excitation/emission of the fluorophore)*

## 4 Analytical Performance

The performance of the FISH Reagent Kit has been established using the procedures described in this instruction for use. Modifications to these procedures may alter the performance of the assay. Any deviations from the instructed protocols must be validated by the user. Users who deviate from recommended test procedures must accept responsibility for interpretation of patient results under these circumstances. Protocol times may vary due to variation in fixation and processing.

## 5 Materials required and supplied with product

| Code #   | Description                                 | Volume     |
|----------|---------------------------------------------|------------|
| LK-095B  | DAPI/Antifade 0.1 µg/ml (DAPI Counterstain) | 0.5 ml     |
| LK-097B  | Antifade (Counterstain diluent)             | 0.5 ml     |
| LK-102B* | 0.4 x SSC / 0.3% Igepal (Wash Buffer I)     | 250 ml     |
| LK-103B* | 2 x SSC / 0.1% Igepal (Wash Buffer II)      | 2 x 250 ml |
| LK-105B  | 2 x SSC/0.5% Igepal                         | 250 ml     |

\*Please note that LK-102 and LK-103 products containing the same number, but different letters are identical except for the volume. A= 100 ml, B=250 ml

## 6 Materials required and not supplied with product

- Ethanol 100%, 85% and 70%
- Fixogum (LK-071A) or rubber cement
- Hot plate with accurate temperature control up to 75 °C or ThermoBrite™ (TS01/TS02)
- Incubator with accurate temperature range of 37 °C
- Water bath with accurate temperature range of 37 - 75 °C
- Plastic and glass coplin jars
- Variable micropipettes (1 µl - 200 µl)
- Fluorescence microscope equipped with suitable filters (see recommendations for Fluorescence Microscopy).

## 7 Warning and precautions

- Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.
- For in vitro diagnostic use. For laboratory professional use only. In case of emergencies check SDS sheets for safety information.
- The kit is intended to use only prior and after to DNA FISH on metaphase and interphase cells from peripheral blood cultures or direct preparations prepared by standard cytogenetic methods.
- For non-automated FISH pre-treatment and hybridization only.
- All materials should be disposed of according to your institution's guidelines for hospital waste disposal.
- It is advised to use the KBI-60006 FISH Digestion kit in case older samples, slides with cytoplasmic background or difficult samples are used or expected.
- For paraffin embedded tissues it is recommended to use pretreatment kits KBI-60004 or KBI-60007.

Labelling according to Regulation (EC) No 1272/2008

| Code #  | Description                                 | Signal Word | Pictogram                                                                           | Hazard Statements |
|---------|---------------------------------------------|-------------|-------------------------------------------------------------------------------------|-------------------|
| LK-095B | DAPI/Antifade 0.1 µg/ml (DAPI Counterstain) | Warning     |  | H315, H319, H335  |
| LK-097B | Antifade (Counterstain diluent)             | Warning     |  | H315, H319, H335  |
| LK-102B | 0.4 x SSC / 0.3% Igepal (Wash Buffer I)     | N.A.        | N.A.                                                                                | N.A.              |
| LK-103B | 2 x SSC / 0.1% Igepal (Wash Buffer II)      | N.A.        | N.A.                                                                                | N.A.              |
| LK-105B | 2 x SSC/0.5% Igepal                         | N.A.        | N.A.                                                                                | N.A.              |

#### Hazard statement(s)

H315 – Causes skin irritation.

H319 – Causes serious eye irritation.

H335 – May cause respiratory irritation

#### Precautionary statement(s) Prevention

P271 – Use only outdoors or in a well-ventilated area.

P261 – Avoid breathing mist/vapours/spray.

P264 – Wash all exposed external body areas thoroughly after handling.

P280 – Wear protective gloves, protective clothing, eye protection and face protection.

#### Precautionary statement(s) Response

P305+P351+P338 – IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P312 – Call a POISON CENTER/doctor/physician/first aider/if you feel unwell.

P337+P313 – If eye irritation persists: Get medical advice/attention.

P302+P352 – IF ON SKIN: Wash with plenty of water and soap.

P304+P340 – IF INHALED: Remove person to fresh air and keep comfortable for breathing.

P332+P313 – If skin irritation occurs: Get medical advice/attention.

P362+P364 – Take off contaminated clothing and wash it before reuse.

Supplemental Hazard Statements: none

Restricted to laboratory professional users

#### Explanation of symbols

|                                                                                     |                                  |
|-------------------------------------------------------------------------------------|----------------------------------|
|   | Contains hazardous substances    |
|  | Country of manufacture           |
|  | Do not use if package is damaged |
|  | Consult electronic IfU (elfU)    |

## 8 Reagent Storage and Handling

Protect DAPI/Antifade from strong and direct light. Store at 2 - 8 °C. The unopened product is stable under these conditions up to the expiry date indicated on the vial label.

Return undiluted reagent to 2 - 8 °C immediately after use.

After opening products are stable for up to 18 Months, but never longer than the expiry date. After this time period the use of the solutions may result in suboptimal results.

## 9 Technical Support

Technical support is available at [www.leicabiosystems.com](http://www.leicabiosystems.com) or via e-mail: [Kreatech.Support@Leicabiosystems.com](mailto:Kreatech.Support@Leicabiosystems.com). Please contact us if desired results are not obtained or any anomalies appear when using this product.

A paper printed copy of the IfU can be requested via e-mail by contacting: [LBSAMS-IFU@leicabiosystems.com](mailto:LBSAMS-IFU@leicabiosystems.com)

The IfU is available in additional languages that can be downloaded from:  
[www.leicabiosystems.com](http://www.leicabiosystems.com)

## 10 Customer Service

Kreatech probes may be ordered through Leica Customer Service +31 20 6919181 or order via e-mail:  
[purchase.orders@leica-microsystems.com](mailto:purchase.orders@leica-microsystems.com).

## 11 History

| IfU Version | Date of issue | Changes                                     |
|-------------|---------------|---------------------------------------------|
| 2.4         | 2016-SEP-27   | Remove the tagline "The pathology Company " |
| 3.0         | 2022-April    | Update according to IVD regulations         |

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