
Tissue Digestion Kit II

REF KBI-60004

Containing sufficient volume for up to 25 or 5x5 tests

Kreatech Biotechnology B.V.
Vlierweg 20
1032 LG Amsterdam
The Netherlands
www.LeicaBiosystems.com

1 Intended use/purpose

Fluorescence in situ hybridization (FISH) identifies or labels target genomic sequences so that their location can be studied. DNA sequences from appropriate chromosome specific probes are first labeled with reporter molecules. The labeled DNA probe is then hybridized to the metaphase chromosomes or interphase nuclei on a slide. After washing, the specimen is screened for the reporter molecules by fluorescence microscopy.

This pretreatment kit is specifically developed to obtain optimal results on heavily cross-linked or difficult formalin-fixed paraffin-embedded tissues (FFPE) (e.g. longer fixation or prolonged storage).

This product is intended for In Vitro Diagnostic Use.

2 Product Description

This pretreatment kit is for non-automated FISH hybridizations specifically developed to use for pretreatment and enzyme digestion of FFPE before DNA FISH hybridization procedure. This kit is specifically developed to obtain optimal results on heavily cross-linked or difficult FFPE tissues. The use of the kit will reverse the cross-links in the tissue, permeabilize the cell membrane and digest cell proteins to ease the binding of the FISH probes with the specific cell DNA.

3 Instructions for use

3.1 Specimens

This pretreatment kit is specifically developed to obtain optimal results on heavily cross-linked or difficult formalin-fixed paraffin-embedded tissues (FFPE) (e.g. longer fixation or prolonged storage).

Specimens' pre-analytics: To obtain a successful FISH the pre-analytics play a major role. The specimen should be correctly fixated and embedded in paraffin wax according to a well-established protocol. Over-fixation or under fixation can cause failure of the FISH procedure or at least low reproducibility/signal intensity. Kreatech™ FISH probes have been tested on FFPE tissue sections which were cut at 3-6 µm. Tissue blocks and cut slides should be stored according to the recommended storage and handling conditions for optimal results, see Leica website. Baking can be done either on the BOND using the *Bake and Dewax protocol or separately in an oven or incubator. Recommended baking temperature and time is up to 16 hrs at 56 °C, or 60 - 120 min at 80 °C.

3.2 Pretreatment

Mount 3 - 6 µm formalin-fixed paraffin-embedded tissue sections on positively charged slides (e.g. aminoalkylsilane). Protocol below is for the use of 5 slides in an upright coplin jar of 50 ml.

1. Bake mounted slides for 2 hours at 80 °C or 16 hours at 56 °C.
2. De-paraffinize warm slides by soaking in xylene or xylene substitute for two times 10 minutes (min).
3. Re-hydrate by soaking in 100%, 85% and 70% ethanol for 3 min each.
4. Wash with dH2O for 3 min at room temperature (RT).
5. Pretreat with 50 ml of 0.2 M HCl (LK-099B) for 20 min at RT, and then wash in dH2O for 3 min at RT.
6. Place slides in 50 ml of Pretreatment Solution B (8% sodium thiocyanate) (LK-100B) at 80 °C for 30 min.
7. Rinse in 50 ml of 2 x SSC (LK-104B) for 3 min at RT.
8. Cover paraffin section with approximately 200 µl Pepsin Solution (LK-101B) and incubate at RT for 5 - 50 min (Time depending on tissue fixation and tissue type. E.g. most breast cancer tissue needs 5 - 15 min digestion; colon tissue 20 - 30 min).
9. Wash in dH2O for 1 min.

10. Wash in 50 ml of 2 x SSC (LK-104B) for 5 min at RT.
11. Dehydrate slides by soaking 70%, 85%, and 100% ethanol for 1 min each time. Air-dry. Proceed with **Probe Preparation**.

3.3 Probe preparation:

Kreatech™ FISH probes for paraffin embedded tissue are supplied Ready to Use (RtU) unless specified otherwise in the product documentation. Consult label on vial and specific probe pack insert for dilution specifics.

3.4 Co-denaturation

Apply 10 µl of probe per 22 x 22 mm field. Cover with glass cover slip and seal with Fixogum or rubber cement. Denature sample and probe on a ThermoBrite (TS-01/02) at 80 ±1 °C for 5 min. Continue with **Hybridization**.

3.5 Hybridization:

Incubate overnight at 37 ±1 °C in a ThermoBrite (TS-01/02) or in a humidified chamber.

3.6 Post-hybridization Wash:

1. Pre-warm 50 ml of Wash Buffer I (0.4 x SSC / 0.3% Igepal) (LK-102) to 72 °C
2. Remove rubber cement.
3. Place slides in 50 ml of Wash Buffer II (2 x SSC / 0.1% Igepal) (LK-103), incubate for 2 min. at RT. Slide off coverslips. Re-use the wash buffer only once for a total of 10 slides.
4. Place slides in 50 ml of pre-warmed Wash Buffer I (0.4 x SSC / 0.3% Igepal) (LK-102), incubate for 2 min at 72 °C (+- 1 °C) without agitation. Re-use the wash buffer only once for a total of 10 slides.
5. Place slides in 50 ml of fresh Wash Buffer II (2 x SSC / 0.1% Igepal) (LK-103), incubate for 1 min. at RT without agitation. Re-use the wash buffer only once for a total of 10 slides.
6. Dehydrate in fresh 70%, 85% and 100% ethanol, incubate for 1 min each at RT. Air dry at RT and proceed to **Counterstaining**.

3.7 Counterstaining:

Apply 15 µl DAPI counterstain (DAPI/Antifade 1 µg/ml) (LK-096B) or a dilution of this using counterstain diluent (LK097B) and apply glass cover slip. Proceed with microscopy.

3.8 Procedural recommendations:

Check protein digestion and pretreatment by applying 15 µl DAPI counterstain and evaluate slides using a fluorescence microscope equipped with a DAPI filter. If the sample is not sufficiently digested: remove cover slip and soak tissue in 2 x SSC (LK-104B) for 2 min and prolong protein digestion. If the sample is over-digested use a fresh sample and reduce protein digestion time.

Temperature and buffer concentration (stringency) of hybridization and washing are important, as lower stringency can result in non-specific binding of the probe to other sequences, and higher stringency can result in a lack of signal. Incomplete denaturation of target DNA can result in lack of signal.

3.9 Recommendations for fluorescence microscopy

For optimal visualization use a well maintained and regularly calibrated microscope equipped with a 100W mercury lamp or other appropriate light source (e.g. Metal Halide Bulb (120W X-cite) or LED based light sources) and a 63x or 100x fluorescence objective. Triple band-pass filters (DAPI/FITC/Cy3 or DAPI/FITC/TRITC) are used to view multiple colors, single band-pass filters are used for individual color visualization. See the Leica website for recommendations on appropriate filters for each microscope brand: <http://www.leicabiosystems.com/ihc-ish/kreatech-fish-probes/>

Suitable excitation and emission range for Kreatech™ fluorophores:

Fluorophore	Excitation	Emission
DAPI	360 ±20 nm	460 ±30 nm
PlatinumBright™ 415	429 ±15 nm	470 ±15 nm
PlatinumBright™ 495	495 ±10 nm	517 ±10 nm
PlatinumBright™ 550	550 ±10 nm	580 ±10 nm

The shown bandwidth corresponds to 80% of the maximum excitation/emission of the fluorophore)

4 Analytical Performance

The performance of Tissue Digestion Kit II has been established using the procedures described in this instruction for use. Modifications to these procedures may alter the performance of the assay. Any deviations from the instructed protocols must be validated by the user. Users who deviate from recommended test procedures must accept responsibility for interpretation of patient results under these circumstances. Protocol times may vary due to variation in fixation and processing. In addition, pepsin incubation time may require optimization depending on processing and fixation conditions.

5 Materials required and supplied with product

Code #	Description	Volume
LK-096B	DAPI/Antifade 1 µg/ml (DAPI Counterstain)	0.5 ml
LK-097B	Antifade (Counterstain diluent)	0.5 ml
LK-099B	0.2 M HCl	250 ml
LK-100B	8% sodium thiocyanate (Pretreatment Solution B)	250 ml
LK-101B	Pepsin Solution	5 ml
LK-102B*	0.4 x SSC / 0.3% Igepal (Wash Buffer I)	250 ml
LK-103B*	2 x SSC / 0.1% Igepal (Wash Buffer II)	2 x 250 ml
LK-104B	2 x SSC	2 x 250 ml

**Please note that LK-102 and LK-103 products containing the same number, but different letters are identical except for the volume. A= 100 ml, B=250 ml*

6 Materials required and not supplied with product

- Xylene
- Ethanol 100%, 85% and 70%
- Fixogum (LK-071A) or rubber cement
- Hot plate with accurate temperature control up to 80 °C or ThermoBrite™ (TS01/TS02)
- Incubator with accurate temperature range of 37 - 56 °C
- Water bath with accurate temperature range of 37 - 80 °C
- Plastic and glass coplin jars
- Variable micropipettes (1 µl - 200 µl)
- Fluorescence microscope equipped with suitable filters (see recommendations for Fluorescence Microscopy).

7 Warning and precautions

- Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.
- For in vitro diagnostic use. For laboratory professional use only. In case of emergencies check SDS sheets for safety information.
- The kit is intended to be use only prior and after to DNA FISH on FFPE tissue. This product is not intended to use for frozen tissue.
- For non-automated FISH pretreatment and hybridization only.
- All materials should be disposed of according to your institution's guidelines for hospital waste disposal.
- The use of the KBI-60007 Tissue Digestion Kit I is recommended for use with conventional FFPE samples.
- For metaphase and interphase cells or other cytological samples it is recommended to use pretreatment kit KBI-60005 or KBI-60006.

Labelling according to Regulation (EC) No 1272/2008

Code #	Description	Signal Word	Pictogram	Hazard Statements
LK-096B	DAPI/Antifade 1 µg/ml (DAPI Counterstain)	Warning		H315, H319, H335
LK-097B	Antifade (Counterstain diluent)	Warning		H315, H319, H335
LK-099B	0.2 M HCl	N.A.	N.A.	N.A.
LK-100B	8% sodium thiocyanate (Pretreatment Solution B)	N.A.	N.A.	N.A.
LK-101B	Pepsin Solution	N.A.	N.A.	N.A.
LK-102B	0.4 x SSC / 0.3% Igepal (Wash Buffer I)	N.A.	N.A.	N.A.
LK-103B	2 x SSC / 0.1% Igepal (Wash Buffer II)	N.A.	N.A.	N.A.
LK-104B	2 x SSC	N.A.	N.A.	N.A.

Hazard statement(s)

H315 – Causes skin irritation.

H319 – Causes serious eye irritation.

H335 – May cause respiratory irritation

Precautionary statement(s) Prevention

P271 – Use only outdoors or in a well-ventilated area.

P261 – Avoid breathing mist/vapours/spray.

P264 – Wash all exposed external body areas thoroughly after handling.

P280 – Wear protective gloves, protective clothing, eye protection and face protection.

Precautionary statement(s) Response

P305+P351+P338 – IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P312 – Call a POISON CENTER/doctor/physician/first aider/if you feel unwell.

P337+P313 – If eye irritation persists: Get medical advice/attention.

P302+P352 – IF ON SKIN: Wash with plenty of water and soap.

P304+P340 – IF INHALED: Remove person to fresh air and keep comfortable for breathing.

P332+P313 – If skin irritation occurs: Get medical advice/attention.

P362+P364 – Take off contaminated clothing and wash it before reuse.

Supplemental Hazard Statements: none

Restricted to laboratory professional users

Explanation of symbols

	Contains hazardous substances
	Country of manufacture
	Do not use if package is damaged
	Consult electronic IfU (elfU)

8 Reagent Storage and Handling

Protect DAPI/Antifade from strong and direct light. Store at 2 - 8 °C. The unopened product is stable under these conditions up to the expiry date indicated on the vial label.

Return undiluted and diluted reagents to 2 - 8 °C immediately after use.

After opening products are stable for up to 18 months, but never longer than the expiry date of each individual product. After this time period the use of the solutions may result in suboptimal results.

9 Technical Support

Technical support is available at www.leicabiosystems.com or via e-mail: Kreatech.Support@Leicabiosystems.com. Please contact us if desired results are not obtained or any anomalies appear when using this product.

A paper printed copy of the IfU can be requested via e-mail by contacting: LBSAMS-IFU@leicabiosystems.com

The elfU is available in additional languages that can be downloaded from: www.leicabiosystems.com

10 Customer Service

Kreatech probes may be ordered through Leica Customer Service +31 20 6919181 or order via e-mail: purchase.orders@leica-microsystems.com.

11 History

IfU Version	Date of issue	Changes
2.4	2016-SEP-27	Remove the tagline "The pathology Company "
3.0	2022 April	Update according to IVD regulations

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