

DDIT3 Break Apart FISH Probe Kit

Introduction

The DDIT3 Break Apart FISH Probe Kit is designed to detect rearrangements in the human DDIT3 gene located on chromosome band 12q13.3. In addition to revealing breaks, which can lead to translocation of parts of the gene, inversion, or its fusion to other genes, the probe set can also be used to identify other DDIT3 aberrations such as deletions or amplifications. Rearrangements and abnormal expression of the DDIT3 gene – also known as CHOP, CEBPZ, CHOP10, CHOP-10, GADD153 – have been observed in myxoid liposarcome and other malignancies.

Intended Use

To detect rearrangements in the human *DDIT3* gene located on chromosome band 12q13.3.


Cont.

Color

LSP DDIT3 5' FISH Probe
LSP DDIT3 3' FISH Probe

CytoGreen
CytoOrange

Probe Design



Not to Scale

LSP DDIT3 5' FISH Probe covers some genomic sequences adjacent to the 5' end of the *DDIT3* gene. LSP DDIT3 3' FISH Probe covers some sequence downstream of the 3' end of the gene. The two probes are flanking sequences across the *DDIT3* gene in which various breakpoints have been observed.

Cat. No.

Volume

CT-PAC057-10-GO

10 Tests (100 µL)

Signal Pattern Interpretation

Normal Patterns

2F*

Abnormal Patterns

Other Patterns

*Overlapping orange and green signals can appear as yellow.

1) Panagopoulos I, et al. *Oncogene*. 12(3):489-94 (1996).
2) Panagopoulos I, et al. *Biochem Biophys Res Commun*. 279(3):838-45 (2000).
3) Pérez-Losada J, et al. *Oncogene*. 19(52):6015-22 (2000).
4) Pérez-Mancera PA, et al. *Oncogene*. 21(11):1679-84 (2002).
5) Rabbitts TH, et al. *Nat Genet*. 4(2):175-80 (1993).

* CE IVD only available in certain countries. All other countries are either ASR or RUO. Please contact your local dealer or our headquarters for more information.