

ENGLISH

For Professional Use Only

MYB-NFIB Fusion/Translocation FISH Probe Kit

Introduction


The MYB-NFIB Fusion/Translocation FISH Probe Kit is designed to detect rearrangements involving the human MYB and NFIB genes, located on chromosome bands 6q23.3 and 9p23, respectively. Rearrangements between the two genes have been observed in adenoid cystic carcinoma and other salivary gland tumors.

Intended Use

To detect rearrangements involving the human *MYB* and *NFIB* genes located on chromosome bands 6q23.3 and 9p23, respectively.

| Cont. | Color |
|------------------------|------------|
| LSP MYB 5' FISH Probe | CytoOrange |
| LSP NFIB 3' FISH Probe | CytoGreen |

Probe Design



LSP MYB 5' FISH Probe covers the entire MYB gene along with upstream (5') and some downstream (3') genomic sequences. LSP NFIB 3' FISH Probe covers the 3' (end) portion of the NFIB gene along with some adjacent genomic sequences. The probe set is optimized to reveal translocations between the two genes.

| Cat. No. | Volume |
|-----------------|-------------------|
| CT-PAC083-10-OG | 10 Tests (100 μL) |

Signal Pattern Interpretation

Normal Patterns **Abnormal Patterns** 202G* Other Patterns

*Overlapping orange and green signals can appear as yellow.

¹⁾ Clappier E, et al. *Blood*. 110(4):1251-61 (2007). 2) Greig KT, et al. *Semin Immunol*. 20(4):247-56 (2008). 3) Ramsay RG & Gonda TJ. *Nat Rev Cancer*. 8(7):523-34 (2008). 4) Murati Å, et al. *Leukemia*. 23(1):85-94 (2009). 5) Di Palma S, et al. *Histopathology*. 64(3):453-5 (2014).

CytoTest Inc. **IVD** 1395 Piccard Drive, Suite 308 -25°C **∕** Rockville, MD 20850, USA

^{*} CE IVD only available in certain countries. All other countries are either ASR or RUO. Please contact your local dealer or our headquarters for more information. DCN032